

Another Innovative Test Designed by WANTAI

WANTAI SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)



(Anterior Nasal Swab & Saliva)

Intended for Self-testing by layperson



WANTAI

Leading infectious diseases
diagnostic test manufacturer

Your Global Partner in IVD



Scan QR code
to watch operation video

Antigen testing has been recognized as a “game changer”, frontline test of the fight against COVID-19. The most commonly used tests are the lateral flow tests which give results in 15min.

An antigen test reveals current coronavirus infection, and it is much faster and convenient to perform than the commonly used RT-PCR test. Anterior nasal and saliva swabs can easily and immediately collect antigens from inside the



Innovative



High Performance




Easy-to-use



Rapid






Sampling:

Anterior Nasal Swab: Insert the absorbent tip of the swab into your nostril. Slowly rotate the swab at least 4 times for a total of 15 seconds. Repeat in your other nostril. **Saliva swabs:** Wipe the upper palate and the inside of the right and left cheeks. Be sure to collect any saliva that may be present on the swab.



Procedure:

- 01 Unscrew the cap of the extraction vial. Place the swab inside, rotate to mix well. **Squeeze the swab** against the inside of the vial to release the liquid from the swab.
- 02 Snap off the end of the swab, leave the swab head in the vial. Screw tightly the cap of the vial and **break the end** of its head.
- 03 Squeeze the vial to add **three (3) drops** of sample to the cassette. Then read the results at **15 minutes**, but no more than 30 minutes.
- 04 Place all the components in the plastic sealable bag, seal it and dispose.
- 05 **Positive results:** Two red lines indicate that SARS-CoV-2 **nucleocapsid antigen** has been detected. **Negative results:** One red line indicates that SARS-CoV-2 antigen has not been detected with this test.



Advanced Diagnostic Technology

- 01 Detection of SARS-CoV-2 Nucleocapsid (N) antigen.
- 02 The sensitivity of the test for all specimen types is 90.11% and the specificity is 99.41%.
- 03 Easy-to-use: works with anterior nasal swab and saliva samples.
- 04 Fast: lateral flow format. Clear answer within only 15 minutes.

Order Information

Components	WJ-2901	WJ-2905	WJ-2910	WJ-2925
Test Cassette	x 1	x 5	x 10	x 25
Extraction Vial	x 1	x 5	x 10	x 25
Disposable Sterile Swab	x 1	x 5	x 10	x 25