




CONFIRMATION TESTING

Wantai confirmation tests are intended for use as supplemental tests for human serum or plasma specimens found to be repeatedly active in HCV and HIV antibody screening procedures.



Recombinant Immunoblot Assays for confirmation of blood-borne pathogens infections

Tests used for the diagnosis of HIV and HCV infection require a high degree of both sensitivity and specificity. This is achieved using an algorithm combining two tests detecting which antibodies. If antibodies are detected by an initial test based on the ELISA method, then a second supplemental test, using the immunoblot procedure confirms or excludes the results of the initial ELISA test. The combination of these two methods is highly accurate.